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
SHIATSU MASSAGER SEAT

PURCHASE ACKNOWLEDGEMENT

JOCCA thanks you for the trust placed in the purchase of our product and we are certain that you will always be satisfied with any of our products.

PRODUCT DESCRIPTION

- a- Massage back rest
- b- Massage Motors
- c- Storing cover for the remote control
- d- Strips to fix the back rest
- e- Massage remote control
- f- Adaptor connector
- g- Adaptor AC/DC

SAFETY WARNINGS

-These instructions are very important. Please read this instructions manual carefully before using the appliance and keep it at hand to refer to it at any time.

-**JOCCA** declines any responsibility for the improper use of the product, or any other use that is not described in this manual.

-Before using the product for the first time, unpack the product and check that it is in perfect conditions. If not, do not use the appliance and contact **JOCCA**, as the guarantee covers any damage from origin or manufacturing defects.

-This product is only indicated for domestic use. Do not use the appliance outdoors. Never leave the appliance exposed to climatic agents such as rain, snow, sun, etc...

-This device is intended for use by adults. Do not allow its use by the persons who do not know how to use it, disabled persons, children or people insensitive to heat.

-Should not be used by children, people with heart problems, thrombosis, serious bodily injury, ultimately if a person has any serious medical problems consult with a doctor.

-Remember that this product has been invented not to treat but to provide relaxation.

-Before connecting the appliance to the mains, check that the cable and plug are in perfect conditions. Also check that the voltage indicated on the appliance is valid for your electrical installation. Plug the deep fryers into a socket with an earthing connection.

-If you detect that the power supply cable or the plug are damaged, do not use the appliance or try to replace the cable yourself. Go directly to the authorised **JOCCA** technical service where they will replace it correctly for you.

-Do not immerse the appliance or any part of it in water or any other liquid.

-To avoid overheating, it is advisable to unroll completely the power cable and to disconnect the plug from the mains when the appliance is not in use.

-If for any reason, the heat that supplies the device is excessive, immediately remove the device from your body.

-Do not pull the cable to unplug the appliance from the mains, always pull from the plug.

-For your safety and for the appliance to work correctly, always use **JOCCA** spare parts and accessories.

-Always unplug a device from electrical network when finish it using, before its cleaning and if it does not work correctly.

-Do not use a device if it is wet; do not leave a device under the exposure of weathering such as rain, sun, snow etc.

-Do not insert pins into a mattress

-Do not use a product if it is wrinkled or folded.

-A prolonged use with maximum level of regulation can lead to burns on the skin.

-Check a device frequently for signs of wear or damage. If such signs exist or a device has been used improperly, consult a distributor before any further use.

-Do not use a device neither with wet or sweat body nor with wet feet or hands.

-When you do not use the device, store it as indicated in the instructions and folding it (the heating seat)

SITUATIONS UNDER WHICH THE BACKREST MUST NOT BE USED

When any of the following situations arises, immediately stop the backrest, switch off the electricity supply, and unplug the appliance:

- When the backrest has been dampened with water, or does not work.
- When the backrest upholstery is damaged or parts of the backrest are exposed.
- When electrical malfunctioning is apparent.
- During a thunderstorm.
- The various parts of the upholstery must be checked to see whether there are any tears or any parts have become unsewn.
- The backrest must not be used if there are any tears in the centre upholstery. The electricity supply must be immediately switched off and unplugged from the mains electricity supply. The use of the massage chair under these circumstances could cause physical damage as a consequence of electric shocks.
- If you feel excessive heat, like a burning back sensation, switch off the backrest immediately. It is recommended you check with your doctor as you may be unaware of suffering from heat sensitivity.

INSTRUCTIONS FOR USE

To switch it on, press the button On / Off  and the pilot light will start flashing and the rollers go down to the lowest position (so they are in the middle of the device now). When they reach the lowest position, the pilot light will stop flashing and remain switched indicating that the device is ready for operation.

Select a function you want to use:

vibration  | rollers  | heating  |

To select the vibration function, press the button Vibration till achieving an intensity level you want. You can select 3 intensity levels of vibration: low, medium and high.

To increase or reduce intensity of vibration, keep pressing the button to alternate from low to medium or high intensity.

For the roller function, you can select an area you want to massage: upper back and neck, lower back or full back, selecting each area by the remote control.

For this first press the button Roller function and then press the button of the area you want. The pilot light will indicate the selected area of massage.

Finally to activate the heating function press the button Heat and the pilot light will be switched on. To switch it off just press the same button.

The device will not be switched off completely until you unplug it from the mains supply.

To switch the device off, press the button On / Off and the pilot light will start flashing, the rollers go down to the lowest position and then up to the middle where they stop and the pilot light will be switched off. If the rollers are in the higher position and we want to switch off the device, the rollers must make the whole way down and half way up until the device turns off.

This procedure is necessary so let the device detecting the exact position of the rollers at each moment.

CLEANING AND STORAGE

-Before cleaning the appliance, make sure it has been disconnected from the mains.

-Do not use abrasive detergents or products that may damage the appliance.

-Clean a device with dry cloth to remove any dust rests or dirt.

-Use dry cloth to clean a remote control.

-Place a device at a clean and dry place.

TECHNICAL SPECIFICATIONS

CHARGER:

110-240V 50/60Hz 1.0A

12V DC 3000mA

This device conforms to the safety requirements and provisions of directives 2006/95/EC on Low tension devices and 2004/108/EC on Electromagnetic Compatibility.

PROTECTION OF THE ENVIROMENT – DIRECTIVE 2002/96/EC

This product has been designed and manufactured with high quality materials and components that can be recycled and reused.

When you see the symbol of a crossed out litter bin on wheels next to a product, this means that the product conforms to the European Directive 2002/96/EC.

Please obtain information about the local recycling system for electric and electronic products.

Follow the local rules and do not dispose of used products by throwing them in the normal bins at home. Proper recycling of your used product will help avoid negative impacts on the environment and people's health.

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